

# Engineers of Confidence.



## Ultrasonic Cleaning The Sound Science

John Nies  
Sr. Product Manager  
1/19/2024

**Belimed**  
Infection Control

# Learning Objectives

Why we use Ultrasonics in the cleaning process?

Understanding the science behind ultrasonics

Knowing what leads to successful ultrasonic cleaning

How to troubleshoot sonic cleaning issues

# What is our cleaning goal?

Remove organic contaminants

- Blood
- Proteins
- Oils
- Tissue

Strip contaminants from all instrument areas

- External & internal surfaces
- Mechanisms

Leave instruments bioburden free, & ready for the next processing step



# Why clean with an Ultrasonic?

Complex instruments are difficult to penetrate

- Box locks
- Drill sleeves
- Flexible reamers
- Graspers
- Cannula
- Robotic tools





# History

**United States Patent** [14] **3,957,252**  
**Storz** [41] **May 18, 1976**

[54] **APPARATUS FOR CLEANING MEDICAL INSTRUMENTS**

[75] **Inventor:** Karl Storz, Tuttlingen, Germany

[73] **Assignee:** Storz-Endoskop GmbH, Schaffhausen, Switzerland

[22] **Filed:** June 17, 1974

[31] **Appl. No.:** 479,822

[36] **Foreign Application Priority Data**

Nov. 7, 1972 Switzerland 1562675  
Apr. 30, 1974 Switzerland 582678

[52] **U.S. Cl.** 259/3 R; 134/1; 134/184; 218/202. 41

[51] **Int. Cl.** B01F 31/02

[58] **Field of Search** 41202; 134/1, 184; 259/1 R, 202. 41; 88/9 R

**References Cited**

**UNITED STATES PATENTS**

2,881,682 3/78 Cooper 4/56  
3,445,992 5/70 Fuchs 1/4/54 R  
3,572,895 3/71 Martin 1/4/53 R

**FOREIGN PATENTS OR APPLICATIONS**


1,002,362 5/55 France 1/4/54

**Primary Examiner—Stanley N. Gilbreth**  
**Assistant Examiner—Alan Cantor**  
**Attorney, Agent, or Firm—Oliver D. Olson**

[57] **ABSTRACT**


Support means is provided for mounting an ultrasonic oscillator for engaging washing water in a conventional sink, for use in cleaning medical instruments. One form of support means is an overflow pipe adapted to be inserted at its lower end in the drain opening of a conventional sink and stabilized at its upper end by an arm extended pivotally from the oscillator power supply control positioned adjacent the sink, the oscillator being mounted on the upper end of the overflow pipe and the overflow pipe having a water inlet opening intermediate its ends for establishing the maximum level of washing water in the sink. Another form of support means is a float structure adapted to float on washing water contained in a conventional sink, the float structure mounting an ultrasonic oscillator for engaging the washing water.

**2 Claims, 2 Drawing Figures**



- Early 1930's - General effect first observed in an RCA test lab
- 1954 – Recognized as a cleaning method by Bendix & others
- 1965 – Becomes the standard for circuit board cleaning (“the RCA clean”)
- Mid-1960's - First medical washing applications

# How The Technology Works



Key component: **PieZoelectric Transducer (PZT)**

- Electronic energy pulses vibrate the PZT's, which emit ultrasound energy
- Ultrasonic energy enters tank water and forms pressure waves



# How The Technology Works



Tiny vacuum voids (bubbles) created


Bubbles implode

Resulting jets create extreme small-scale effects:

- 20,000 PSI pressure
- 5,000K temperature

Effect is called Cavitation

# What does “Cavitation” do?



Creates a powerful effect that blasts off heavy soil

- Shear forces, as bubbles collapse, rip particulates off surfaces
- Ultra = beyond, Sonic = sound. Sound beyond what you can hear
- High frequency (>40khz) = smaller bubbles & fine detail cleaning
- Electronics parts, hard disk drives, optics
- Low frequency(40khz) = larger bubbles & more “scrubbing” power
- Medical instruments, automotive parts

Effect isn't vibration-based, it's all about the cavitation



# Where else have we seen Sonics?




# Other Applications

Penetrates crevices

- Anywhere liquid can reach
- Areas inaccessible by pressure alone

Many uses

- Automotive parts
- Electroplating
- Pharmaceutical production
- Electronics components
- Optics
- And (of course), Healthcare



# To Sonic...or not to Sonic

Which instruments can I clean this way?

- Consult the instrument IFU
- It will indicate if ultrasonic cleaning is needed

Which instruments are best served by sonication?

- The more complex, the better the “fit”
- The less complex, the less added value sonics provide





# Avoid the Technology Limitations

Only works in liquid

- Always fully submerge instruments

Is impeded by solid masses

- Avoid solid metal trays
- Don't use solid/perforated plastic trays
- Best tray type – open-weave wire mesh or heavily perforated metal
- Don't stack trays



# Power Concentration / Density

Power density: transducer power / tank water volume

- Example: 1000W power / 15 gallon tank = 66W/Gal.
- Higher the ratio, the more ultrasonic power is being focused within the tank
- Minimum acceptable is generally defined as 48-50 W/Gal.
- Power output typically fixed, can't "turn it up to 11"
- Sonics with multiple fill levels may have multiple power generators

**Recommendation:** use only enough water to fully cover instruments



# Unit loading



Should you fill the chamber up to the top with instruments?

- Most manufacturers specify a tray capacity which is lower than the max water level
- Overloading often extends cleaning duration
- Overloading creates risk – will this load be clean? Will all instruments be submerged?
- Know the lift weight capacity of the sonic



# Water Quality

AAMI TIR34 lists “key water quality factors to consider for ultrasonic cleaning”:

- a) Water hardness
- b) Water temperature
- c) Ionic contaminants (eg. Chlorides & heavy metals)
- d) Microbial level
- e) Bacterial endotoxins

Two areas of concern

- Unit function (scale deposits)
- Wash quality (all the rest)



# Chemistry



Use non-foaming, sonic-indicated, detergent

Maintain water within detergent manufacturer's specified temperature range

Dose correctly

- Concentration
- Some sonics dose by time, some by flowmeter
- Flowmeter dosing is more accurate

# Wash Time



There is no one “standard cycle”

Most cycles specified by instrument manufacturer’s IFU

IFU cycles can range from 8 to 20 minutes

IFU cycles seem to assume minimum loading – heavy cycles may need more time



# Degassing



Sonic cavitation much less effective when air is dissolved in water

“Degassing” is forcing the dissolved air to escape the water in the tank

Can be manual (add degas time to cycle time), or automatically programmed into machine

Longer degas when tank is freshly filled

Shorter for additional uses of the same water

# Tank Water Quality



How often should tank water be changed?

**AAMI ST79:2017, Section 7.6.4.4.1 d)**

**“...after each use (a “use” should be defined in the health care facility’s policies and procedures):”**



# Lumen flushing



Pushes water through lumens to remove released bioburden

Prevents efficiency-robbing air pockets within lumen

If you're washing cannulated instruments, you need to have this feature (DaVinci)

Especially effective when used both in wash and rinse cycles





# Rinsing



## Debris removal – surfaces

- Post-sonic rinsing useful to prevent debris re-deposition & soap build-up

## Do you need to rinse?

- Are you bypassing the washer/ disinfector?
- Are you processing lumened instruments?
- Do you have IFU's that specify post-sonic rinse?
- Processing daVinci Endowrists?

If yes to any of these, you may need a sonic with rinse

# Purified (DI/RO) Water Rinsing

Do you need to DI rinse?

- Is it specified in the instrument IFU?
- Are your instruments bypassing the washer/disinfector?
- If No to both & the instruments don't have a lumen, then probably not.

AAMI TIR-34 2007 implies treated water is a “factor to consider” for automated cleaning sonics, but doesn't explicitly spell this out



# Verifying Sonic Function

Is my sonic creating cavitation (is it working?)

“Tin Foil” test


- Traditional method, inexpensive
- Proves that transducers are functioning

Pre-packaged vial

- Coated beads that change color in response to sonic energy
- More expensive, yet easier to interpret

AAMI ST79 considers a sonic a “mechanical washer”

- Test daily





# Verifying Sonic Cleaning

Is my functioning sonic having the desired effect?

Test coupons

- Backing material, coated with challenge media
- Some use blood & fibrin, others ink dots or other material
- Very similar to washer test coupons
- Make sure that all test surfaces are immersed when in use
- Check with sonic manufacturer for tips/recommendations before using



# Resolving Wash Issues

“My instruments aren’t exiting the sonic clean”, or “I’m failing sonic validation testing”

## Check:

- Chemistry concentration & type (tune the power generators)
- Tank water temperature
- Water quality, especially hardness
- Degas time
- Unit overloading
- Lumen flush connected
- Instrument rinse
- Basic sonic function

Adjust / change only one of the above at a time to pinpoint the true cause.

Run a few test loads to double-check cycle settings



# Recording Cycle Results

Still uncommon for sonics

Cycle parameter output types

Paper printout

- Traditional approach, found in some older machines
- Difficult interface with electronic record-keeping.

USB

- Easy commercially available storage
- Stand-alone, no automatic link with network resources

Networked

- Seamless connection to instrument-tracking systems





# Instrument Drying



Some older sonics offer a dry feature

- Post-sonic drying useful for direct-to-low-temp instruments

Do you need to dry?

- Do you/will you sonic low-temp (H<sub>2</sub>O<sub>2</sub>/Ozone) instruments?
- Will those low-temp instruments bypass the washer/disinfector?
- Consider a drying cabinet – more efficient

# Disinfection

Some sonics offer a thermal disinfection feature

Do you need to disinfect within the sonic?

- Are your instruments bypassing the washer/disinfector?




# Sonics and daVinci

STEP 14 (optional)

## Thermal Disinfection

**A**



**DISINFECT**  
Perform disinfection per hospital policy and regional guidelines. Thermal disinfection is not a substitute for reprocessing.

| PARAMETER   | VALUE                           |
|-------------|---------------------------------|
| Temperature | 185 °F - 199 °F (85 °C - 93 °C) |
| Time        | 1 - 5 minutes                   |

Thermal disinfection is OPTIONAL


Always get latest IFU and for Endowrists from Intuitive's Sterile Reprocessing Support Team (SRS) Rep:

To schedule an In-Service at your hospital, contact  
**Customer Service** at **1.800.876.1310** or e-mail **SRS-Support@intusurg.com**  
 please provide your Hospital Name, City/State & Contact Information

- Endowrists can be run through a thermal disinfect cycle in the washer, in a basket. ALL previous washing steps MUST be completed per the IFU.
- Time and temp of WD Thermal Disinfect step can be set within the IFU ranges.



# Ergonomics




Look for a tray/lid lift function

- Reduces the strain of lifting heavy instrument sets out of a deep chamber

Should I worry about ergonomics & sonics?

- Do your sonics have a chamber that's more than 10" deep (measured from the top work-surface)?
- Are you sonic'ing full, heavy sets?
- If yes to either, then a powered elevator will be a good choice
- Check loading height for ease of use

# Balance Workflow between Sinks and Sonic



Where should sonics be in my SPD?

- As close to the processing sink as possible
- If consistently sonic'ing, long path to the sonic can cost time, and interfere with other SPD staff
- If sonic'ing rarely, walking over to a central sonic is more acceptable

# General Upkeep

Keep lid closed when not in use to limit evaporation

Inspect the detergent reservoir

Clean tank and exterior panels at least once a week

Inspect plumbing & hoses for leaks

Any additional maintenance as called out in the IFU





# Major Maintenance Watchouts

## Transducers

- Constant vibration can cause failure in epoxy bonds, brazing is more reliable
- Tank bond failure results in loss of function


## Mechanicals

- Lid/tray lift mechanisms can wear out over time
- The more complex the mechanism, the more likely it is to fail

## Electronics

- Keep water away from displays & controls, particularly if below tank opening

Service support is key



# Summary



Sonics are a powerful and appropriate soil removal tool



Sonic setup is key to good results



Reliable tools exist to validate and verify unit function



There's no "one-size fits all"; select a sonic based on your expected usage

# Questions?

Thank-You for your time today.



**Beli**med  
Infection Control



# Contact

Belimed Inc. is a leading supplier of medical and surgical instrument sterilization, disinfection, and cleaning products and services.

With 10 Belimed companies throughout Europe, North America and China, Belimed is represented worldwide by a strong network in over 80 countries.

## USA Headquarters

Belimed, Inc.  
8351 Palmetto Parkway, Suite 101  
Ladson, SC 29456  
1.800.451.4118  
[info.us@belimed.com](mailto:info.us@belimed.com)

## Global Headquarters

Belimed AG  
Grienbachstrasse 11  
6300 Zug  
Switzerland  
+41 41 449 78 88  
[info@belimed.com](mailto:info@belimed.com)

REV1.2022 - © 2022 Belimed, Inc.



Engineers of Confidence.

**Belimed**  
Infection Control

For more information  
visit [www.belimed.com](http://www.belimed.com)

