STERIS UNIVERSITY



Keeping Biofilms At Bay In Your Reprocessing Loop

Disclosures

- Successful completion: Participants must complete the entire program and submit required documentation. No partial credit will be given.
- 2. Conflict of interest: Employee of STERIS.
- 3. Commercial company support: Fees are underwritten by education funding provided by STERIS.
- 4. Non-commercial company support: None.
- 5. Alternative/Complementary therapy: None.

Continuing Education

- STERIS Corporation is an approved provider of continuing nursing education by CBRN – provider # CEP 11681 and an approved Administrator Education Unit (AEU) and Infection Prevention Control (IPCH) provider by BASC – provider # 1417.
- This program is approved for:
- O hour(s) of GI Specific content credit by ABCGN (American Board of Certification for Gastroenterology Nurses),
- 1.5 AEU(s) & 1.5 IPCH(s) by BASC (Board of Ambulatory Surgery Certification), and
- 1.5 contact hour(s) of continuing education credit
- CBRN (California Board of Registered Nursing);
- CBSPD (Certified Board for Sterile Processing and Distribution); and
- HSPA (Healthcare Sterile Processing Association).

Continuing Education

- STERIS Corporation is an approved provider of continuing nursing education by CBRN
 – provider # CEP 11681 and an approved Administrator Education Unit (AEU) and
 Infection Prevention Control (IPCH) provider by BASC provider # 1417.
- · This program is approved for:
- <u>0</u> hour(s) of GI Specific content credit by ABCGN (American Board of Certification for Gastroenterology Nurses),
- 1.5 AEU(s) & 1.5 IPCH(s) by BASC (Board of Ambulatory Surgery Certification), and
- <u>1.5</u> contact hour(s) of continuing education credit
- CBRN (California Board of Registered Nursing);
- CBSPD (Certified Board for Sterile Processing and Distribution); and
- HSPA (Healthcare Sterile Processing Association).



Learning Objectives

- List three requirements for biofilm formation
- · Describe the hazards of medical biofilm formation
- Identify areas where biofilms hide and ways to mitigate their formation













Mature Biofilms Are Tough

- Thick protective matrix
- Increased Resistance to:
- Cleaning Chemistries
- Disinfectants
- Sterilants
- Antibiotics







Part of Quality Management System (QMS) Part of Quality Management System (QMS) Risk Management Plan Identify Mitigate Monitor

















Mitigation - Water

- Silver and Copper Treatment
- Remove pipe dead legs
- Regular decontamination of holding tanks
- Disinfect supply lines
- Water Treatment Filter Preventative Ma

<u>CMS</u>

"water management Policies and procedures to reduce the risk of growth and spread of Legionella and other opportunistic pathogens in building water systems."











| Point Of Use Treatment | | |
|-------------------------|--------------|--------------------------|
| Free of Gross Debris | Flush Lumens | Prepare for Transport |
| | | |



Mitigation – Point-Of-Use And Transport

- Increase transport and cleaning capacity
- Identify and correct pretreatment misses
- Staff Education & competency
- Product availability
- Transport device cleaning procedures













Mitigation – Manual Cleaning

- Monitor soiled device wait times
- Follow device manufacturer's written IFU
- Delayed cleaning protocol
- Cleaning verification tests
- Lighted magnification
- Borescope





<section-header><section-header>





Mitigation – Automated Cleaning

- Clean, decontaminate and disinfect equipment
- Inspect, decontaminate and disinfect drains
- Repair Leaks

Cleaning Chemistries





Mitigation – Cleaning Chemistries And Accessories

- Temperature and humidity control
- No Topping Off
- Clean & decontaminate accessories
- Consider single use items
- Clean, dry, and HLD or sterilize water bottles









Methods Of Manual Drying

- Instrument air purging
- Alcohol flush
- Non-linting towels



Mitigation - Drying

- Manufacturer's written IFU
- Pass through
- Drying cabinets / compressed air
- Inspection: Lighted Magnifiers / Borescopes / Moisture Check



Assembly/Prep And Pack Area Moisture & Porous Materials Linen & Towels

Assembly / Prep And Pack - Mitigation

- Clean & decontaminate
- Control humidity levels
- · Replace damp towels & pads
- No lotions, food and drinks
- Inspect











Mitigation – Steam Sterilization

- Follow device manufacturer's written IFU
- Cool loads completely
- Steam quality testing
- Preventative maintenance
- Inspect
- Clean & decontaminate floor drains



Mitigation - Liquid Chemical Sterilization

- Clean & decontamination
- Clean up drips
- Perform regular maintenance per the manufacturer's written IFU
- Inspect









High-Level Disinfection - Mitigation

- Follow manufacturer's written IFU
- Clean & disinfect AER and accessories
- Preventative maintenance
- Expiration
- Inspect
- Concentration testing







Mitigation - Storage

- Ensure dry devices
- No storage beneath water lines
- Control humidity
- Clean & disinfect endoscope cabinets
- Change pads
- Inspect for moisture

















Mitigation – Loaned Instrument Trays

- · Timely receipt
- Follow Manufacturer's Written IFUs
- Inspect
- Follow Manufacturer's Written Dry times
- Cool trays

Summary

- · Three requirements for biofilm formation
- Describe the hazards of medical biofilm formation
- · Identify areas where biofilms hide and ways to mitigate their formation

Action Plan

- · Become familiar with your SPD biofilm risks
- Update risk management plan
- Develop audit checklists
- Establish audit schedule

References

- Alfa, MJ. (Spring, 2012) The "Pandora's Box" Dilemma: Reprocessing of Implantable Screws and Plates in Orthopedic Tray Sets. The Society for Healthcare Epidemiology of America Horizons. pp 55-59.
- Association for the Advancement of Medical Instrumentation. (2018). ANSI/AAMI ST58:2013/(R)2018 Chemical sterilization and high-level disinfection in health care facilities. Arlington, VA: Author.
- Association for the Advancement of Medical Instrumentation. (2020). ANSI/AAMI ST79:2017 & 2020 Amendments A1, A2, A3, A4 (Consolidated Text) Comprehensive guide to steam sterilization and sterility assurance in health care facilities. Arlington, VA: Author.
- Association for the Advancement of Medical Instrumentation. (2017). ANSI/AAMI ST90:2017 Processing of heath care products—Quality management systems for processing in Health care facilities. Arlington, VA: Author.

References

- Association for the Advancement of Medical Instrumentation. (2014/(R) 2017. ANSI/AAMI TIR34: (2014/(R) 2017 Water quality for reprocessing of medical devices. Arlington, VA: Author.
- Association of periOperative Registered Nurses. (2020). Guidelines for perioperative practice. Denver, CO: Author.
- Donlan, Rodney. (2002) Biofilms: Microbial Life on Surfaces. Emerging Infectious Diseases 8:881-890. Centers for Disease Control and Prevention
- Hoiby N. A short history of microbial biofilms and biofilm infections. APMIS. 2017 Apr;125(4):272-275. doi: 10.1111/apm.12686. PMID: 28407426.A

Questions?

Obtain your Certificate

- Thank you for attending this CE activity
- Please complete and submit the evaluation form

Thank you!

Visit us at university.STERIS.com

Supporting a Strong Department Developing Future Leaders

UNIVERSITY